
NMRViewer 1.0 - (C) INRA UMR 1332 Daniel Jacob

NMR spectra viewer: An overview of how it works




Web design based on the jQuery framework (client side), PHP, C, gnuplot 4.3 (server side)

NMR spectra viewer: Why it runs so fast?

Because **Text to Binary conversion is time consuming**, so we do not manipulate data as text format but **only as binary format in foreground process**.

Text to Binary conversion is set up in a background process, only once.



Trade off between consistency and efficiency

NMR spectra viewer: How to process with an ISA-TAB archive?

Step 1: unzip the isatab archive

```
[djacob@proticdb nmrview]$ 11 ./data/HUGO07/isatab/
-rw-r--r- 1 diacob users 317389 Mar 26 14:56 BB BBI 08M07 melon07 13091.zip
-rw-r--r 1 djacob users 317361 Mar 26 14:56 BB BBI 08M07 melon07 13092.zip
-rw-r--r- 1 djacob users 318559 Mar 26 14:56 BB BBI 08M07 melon07 13101.zip
-rw-r--r- 1 diacob users 318787 Mar 26 14:56 BB BBI 08M07 melon07 13102.zip
-rw-r--r- 1 djacob users 318251 Mar 26 14:56 BB BBI 08M07 melon07 13111.zip
-rw-r--r- 1 djacob users 318122 Mar 26 14:56 BB BBI 08M07 melon07 13112.zip
-rw-r--r- 1 diacob users 309234 Mar 26 14:56 BB BBI 08M07 melon07 13121.zip
-rw-r--r- 1 djacob users 318059 Mar 26 14:56 BB BBI 08M07 melon07 13122.zip
-rw-r--r- 1 djacob users 318926 Mar 26 14:56 BB BBI 08M07 melon07 13131.zip
-rw-r--r- 1 djacob users 309249 Mar 26 14:56 BB BBI 08M07 melon07 13132.zip
-rw-r--r-- 1 djacob users 317804 Mar 26 14:56 BB BBI 08M07 melon07 13141.zip
-rw-r--r-- 1 djacob users 318394 Mar 26 14:56 BB BBI 08M07 melon07 13142.zip
-rw-r--r- 1 djacob users 323837 Mar 26 14:56 BB BBI 08M07 melon07 13151.zip
-rw-r--r-- 1 djacob users 325276 Mar 26 14:56 BB BBI 08M07 melon07 13152.zip
-rw-r--r-- 1 djacob users 324708 Mar 26 14:56 BB BBI 08M07 melon07 13161.zip
-rw-r--r- 1 diacob users 323442 Mar 26 14:56 BB BBI 08M07 melon07 13162.zip
-rw-r--r- 1 djacob users 324762 Mar 26 14:56 BB BBI 08M07 melon07 13171.zip
-rw-r--r- 1 djacob users 323527 Mar 26 14:56 BB BBI 08M07 melon07 13172.zip
-rw-r--r- 1 djacob users 324286 Mar 26 14:56 BB BBI 08M07 melon07 13181.zip
-rw-r--r- 1 djacob users 322688 Mar 26 14:56 BB BBI 08M07 melon07 13182.zip
-rw-r--r- 1 djacob users 323630 Mar 26 14:56 BB BBI 08M07 melon07 13191.zip
-rw-r--r- 1 djacob users 324020 Mar 26 14:56 BB BBI 08M07 melon07 13192.zip
-rw-r--r-- 1 djacob users 324789 Mar 26 14:56 BB BBI 08M07 melon07 13201.zip
-rw-r--r- 1 djacob users 324475 Mar 26 14:56 BB BBI 08M07 melon07 13202.zip
-rw-r--r- 1 diacob users 311251 Mar 26 14:56 BB BBI 08M07 melon07 13331.zip
-rw-r--r- 1 djacob users 316116 Mar 26 14:56 BB BBI 08M07 melon07 13332.zip
-rw-r--r 1 djacob users 311180 Mar 26 14:56 BB BBI 08M07 melon07 13341.zip
-rw-r--r- 1 djacob users 314663 Mar 26 14:56 BB BBI 08M07 melon07 13342.zip
-rw-r--r- 1 djacob users 319548 Mar 26 14:56 BB BBI 08M07 melon07 13351.zip
-rw-r--r- 1 djacob users 317141 Mar 26 14:56 BB BBI 08M07 melon07 13352.zip
-rw-r--r-- 1 djacob users 326022 Mar 26 14:56 BB BBI 08M07 melon07 13361.zip
-rw-r--r- 1 djacob users 324277 Mar 26 14:56 BB BBI 08M07 melon07 13362.zip
-rw-r--r- 1 djacob users 320622 Mar 26 14:56 BB BBI 08M07 melon07 13371.zip
-rw-r--r- 1 djacob users 310643 Mar 26 14:56 BB BBI 08M07 melon07 13372.zip
-rw-r--r- 1 djacob users 319364 Mar 26 14:56 BB BBI 08M07 melon07 13381.zip
-rw-r--r- 1 djacob users 319049 Mar 26 14:56 BB BBI 08M07 melon07 13382.zip
-rw-r--r- 1 djacob users 8414 Mar 26 14:56 a hugo 2007 metabolite profiling NMR spectroscopy.txt
-rw-r--r 1 djacob users 10906 Mar 26 14:56 i Investigation.txt
-rw-r--r-- 1 djacob users 5023 Mar 26 14:56 s Hugo 2007.txt
```

Step 2: Preparing for viewing

BB BBI 08M07 melon07 13382 (36/36) ...

NAME = 07.HU3.6.1338.2

Pack ...OK

```
[djacob@proticdb nmrview]$ sh ./scripts/buildcsv.sh -d ./data/HUG007/isatab/ -o samples.csv -v
i Investigation.txt ... OK
s Hugo 2007.txt ... OK
a hugo 2007 metabolite profiling NMR spectroscopy.txt ... OK
Extract columns from Study Assav File : 1,28 ... OK
Factors: ; Harvest date
Extract columns from Study File: 9,10 ... OK
Merging ...OK
[djacob@proticdb nmrview] $ sh ./scripts/init bin.sh -i ./data/HUG007/isatab/samples.csv -o ./data/HUG007/procs/ -h
BB BBI 08M07 melon07 13091 (1/36) ...
NAME = 07.HU1.1.1309.1
BB BBI 08M07 melon07 13092 (2/36) ...
NAME = 07.HU1.1.1309.2
BB BBI 08M07 melon07 13101 (3/36) ...
NAME = 07.HU1.2.1310.1
_____
                                       [djacob@proticdb nmrview]$ 11 ./data/HUGO07/procs/
BB BBI 08M07 melon07 13102 (4/36) ...
                                       total 9756
NAME = 07.HU1.2.1310.2
                                        -rw-r--r-- 1 djacob users 28 Mar 26 15:03 factors
OK
                                       -rw-r--r-- 1 djacob users 1812 Mar 26 15:03 samples.csv
                                        -rw-r--r-- 1 djacob users 9961496 Mar 26 15:04 specs.pack
BB BBI 08M07 melon07 13372 (34/36) ...
NAME = 07.HU3.5.1337.2
BB BBI 08M07 melon07 13381 (35/36) ...
NAME = 07.HU3.6.1338.1
```

Step 3: configuration file : ./etc/nmrproc.conf

```
# Declaration of global variables
# Warning: Take care not putting spaces each side of equal sign, in order to be executed as a script shell
# VARNAME=anything including spaces

DEBUG=0
# The root of the bulk of your installation.
ROOTDIR=/mnt/var/sys/projets/nmrview
# The root of the datasets
DATASETS=/mnt/var/sys/projets/nmrview/data

BASEURL=www.cbib.u-bordeaux2.fr
URLROOT=/SPECNMR
```

jQuery UI plugin: simple.html

```
<!DOCTYPE html PUBLIC "-//W3C//DTD XHTML 1.0 Strict//EN" "http://www.w3.org/TR/xhtml1/DTD/xhtml1-strict.dtd">
<html xmlns="http://www.w3.org/1999/xhtml" xml:lang="en" lang="en">
<head>
  <meta charset="utf-8" />
 <TITLE>NMRProcFlow - Viewer</TITLE>
  <link rel="stylesheet" href="/SPECNMR/css/redmond/jquery-ui.min.css">
  <script src="http://code.jquerv.com/jquerv-1.9.1.js"></script>
  <script src="http://code.jquerv.com/ui/1.10.0/jquerv-ui.js"></script>
  <script src="/SPECNMR/js/ui.nmrview.js"></script>
 <script>
  $(function() {
     $( "#nmrview" ).nmrview({
       url: '/SPECNMR/view', dataset: 'HUGO07', colorby: 0, width: 'auto'
     }).show();
 });
  </script>
</head>
<body>
  <div id="nmrview"></div>
</body>
</html>
```

